



PROVOKE®

RAT MONITORING GEL

**Non-Toxic • Non-Staining
Hypoallergenic**



NET WT:
8 oz (227 g)

- PROVOKE® can be used with any type of monitoring device or mechanical trap... plastic, wood, live catch, etc.
- Remove top and apply a small amount of PROVOKE® onto appropriate area of device or trap.

Manufactured by: Bell Laboratories, Inc.
3699 Kinsman Blvd., Madison, WI 53704 U.S.A.
www.belllabs.com | **MADE IN USA**



102621
523PR-3